

9025 Technology Dr. • Fishers, IN 46038-2886
 800.387.0672 • 317.570.7020 • Fax 317.570.7034
 info@bangslabs.com • www.bangslabs.com



B E A D S ● A B O V E T H E R E S T™

DESCRIPTION

The BioMag[®] Multi-6 Microcentrifuge Tube Separator is a durable, plastic magnetic separation unit designed to accommodate six (6) 1.5mL microcentrifuge tubes. The magnets and tube holders are positioned to efficiently magnetically separate microliter volumes common in molecular biology applications.

PROCEDURE

Magnetic separations can be performed on six (6) 1.5mL microcentrifuge tubes simultaneously. Once the tubes are placed in the tube holders, the magnetic particles are pulled to the side of the unit to complete a separation in less than 1 minute. Upon completion of the separation, the supernatant may be aspirated, leaving the magnetic particles as a pellet on the side of the tube.

STORAGE AND STABILITY

The BioMag[®] Multi-6 Microcentrifuge Tube Separator should be wiped clean with mild detergent. Do not immerse. Harsh chemicals and organic solvents may damage the surface of the unit.

SAFETY

Strong magnetic field. This separation unit should be used and stored a safe distance from magnetic media, such as computers, videotapes, credit cards, etc.

This product is for research use only and is not intended for use in humans or for *in vitro* diagnostic use.

ORDERING INFORMATION

Cat. Code	Description	Size
MS002	BioMag [®] Multi-6 Microcentrifuge Tube Separator	1 each

Order online anytime at www.bangslabs.com.

