


Colloidal Arrays & Nanosphere Lithography



A variety of applications exist in which uniform spherical nano- or microparticles are deposited onto a solid substrate or in some type of 3-D matrix in a highly ordered fashion, e.g. as self-assembled monolayers or close-packed lattices for creating equidistant spacing^{1,2}, nanosphere lithography^{3-7,19}, colloidal crystal assembly,⁷⁻¹³ generating photonic crystals¹⁴⁻¹⁶, fabricating gradient surfaces,¹⁷ biomolecular micropatterning,¹⁸ etc. After bead deposition some of these applications involve the filling or modification of interstices with a material that solidifies. The beads are then fully removed or otherwise altered in order to create highly uniform pores or patterned nanostructures in the solidified material. These modified materials find utility in various optical, semiconductor, and biosensing applications.^{1-17,20}

Numerous techniques exist for depositing beads such as: 2D self-assembly at double- or triple-phase boundaries, spin-coating, ink-jet printing, sedimentation (via gravity, centrifugation, or filtration), convective assembly and electrophoretic methods, two-substrate vertical deposition via confinement between thin films, slit filling, and physical methods such as compression or compaction of dried particles, etc.¹⁻¹⁸ Likewise, a variety of methods have been employed to fill interstices (e.g. polymerization or sol-gel hydrolysis of a liquid precursor, electroless or electrochemical deposition, precipitation, deposition of smaller colloidal particles⁷) or to modify mask openings in nanosphere lithography (e.g. thin-film deposition, sputter erosion, and ion implantation⁴). Techniques for partially or fully removing beads in order to create uniform pores include oxygen plasma treatment, calcination, and chemical dissolution, e.g. organic solvents for polymer beads, strong acid for silica beads, etc.^{5,8,10}

Due to their stringent size and shape uniformity, our polymer, and silica microspheres have been widely used in the types of applications cited above.^{3,11,13,15,16} Ranges for standard bead diameters offered and our complete inventory of available lots, as well as a wide variety of additional microsphere products may be found at www.bangslabs.com.

PS, PS COOH: ~25nm – 90µm
silica: ~150nm – 5µm

REFERENCES:

- Buxboim A, Rajagopal K, Brown AEX, Discher DE. (2010) How deeply cells feel: methods for thin gels. *J Phys Condens Matter*. May 19;22(19):194116.
- Monobe H, Awazu K, Shimizu Y. (2006) Alignment Control of a Columnar Liquid Crystal for a Uniformly Homeotropic Domain Using Circularly Polarized Infrared Irradiation. *Advanced Materials*. March;18(5):607-10.
- Kuo C-W, Shiu J-Y, Chen P. (2003) Size- and Shape-Controlled Fabrication of Large-Area Periodic Nanopillar Arrays. *Chem. Mater.*;15(15):2917-20.
- Lindner J. (2013) Nanosphere lithography for device fabrication. *SPIE Newsroom*. Oct. 16.
- Hua F, Lvov Y, Cui T. (2002) Spatial patterning of colloidal nanoparticle - based thin film by a combinative technique of layer-by-layer self-assembly and lithography. *J Nanosci Nanotechnol*. Jun-Aug;2(3-4):357-61.
- Yan Q, Chen A, Chua SJ, Zhao XS. (2006) Nanosphere lithography from template-directed colloidal sphere assemblies. *J Nanosci Nanotechnol*. Jun;6(6):1815-8. PMID: 17025090
- Zhang J, Li Y, Zhang X, Yang B. (2010) Colloidal Self-Assembly Meets Nanofabrication: From Two-Dimensional Colloidal Crystals to Nanostructure Arrays. *Advanced Materials*;22(38):4249-4269.
- Velev OD, Lenhoff AM. (2000) Colloidal Crystals as templates for porous materials. *Current Opinion in Colloid & Interface Sciences*;5:56-63.
- Mihi A, Ocaña M, Míguez H. (2006) Oriented Colloidal-Crystal Thin Films By Spin-Coating Microspheres Dispersed in Volatile Media. *Advanced Materials*;18(17):2244-49.
- Sun Z, Li Y, Wang Y, Chen X, Zhang, J, Zhang K, Wang Z, Bao C, Zeng J, Zhao B, Yang B. (2007) Three-dimensional colloidal crystal-assisted lithography for two-dimensional patterned arrays. *Langmuir*. Oct 9;23(21):10725-31.
- Lumsdon SO, Kaler EW, Velev OD. (2004) Two-Dimensional Crystallization of Microspheres by a Coplanar AC Electric Field. *Langmuir*. 20;2108-16.
- Sun J, Tang C-j, Zhan P, Han Z-l, Cao Z-S, Wang Z-L. (2010) Fabrication of Centimeter-Sized Single-Domain Two-Dimensional Colloidal Crystals in a Wedge-Shaped Cell under Capillary Forces. *Langmuir*. ;26(11):7859-64.
- King SB, Dorfman KD. (2013) Role of Order during Ogston Sieving of DNA in Colloidal Crystals. *Anal. Chem.*;85(16):7769-76.
- Yamasaki T, Tsutsui T. (1998) Spontaneous emission from fluorescent molecules embedded in photonic crystals consisting of polystyrene microspheres. *Applied Physics Letters*. April 20;72(16):1957-9.
- Inam FA, Gaebel T, Bradac C, Stewart L, Withford MJ, Dawes JM, Rabeau JR, Steel MJ. (2011) Enhanced spontaneous emission from nanodiamond colour centres on opal photonic crystal. *New J. Phys.*;13:073012.
- Park J, Moon J, Shin H, Wang D, Park M. (2006) Direct-write fabrication of colloidal photonic crystal microarrays by ink-jet printing. *Journal of Colloid and Interface Science.*;298:713-9.
- Zhang J, Xue L, Han Y. (2005) Fabrication Gradient Surfaces by Changing Polystyrene Microsphere Topography. *Langmuir*;21:5-8.
- Wang C, Yap FL, Zhang Y. (2005) Micropatterning of polystyrene nanoparticles and its bioapplications. *Colloids Surf B Biointerfaces*. Dec 30;46(4):255-60. PMID: 16356699
- Burmeister, F., Schäfle, C., Matthes, T., Böhmisch, M., Boneberg, J., & Leiderer, P. (1997). Colloid monolayers as versatile lithographic masks. *Langmuir*, 13(11), 2983-2987.
- Velev, O. D., & Kaler, E. W. (1999). In situ assembly of colloidal particles into miniaturized biosensors. *Langmuir*, 15(11), 3693-3698.