


Europium Chelate

Stable, high-intensity Eu microspheres for time-resolved and other fluorescence assays.

Lateral flow and other rapid tests fulfill important roles in the diagnostic landscape. Affordable and easy to use, they are particularly important for delivering diagnostic capability to programs with critical need, limited resources, or remote / de-centralized laboratories. While many significant tests have been developed using conventional particles (e.g. colloidal gold), the use of europium chelate (Eu(III)) nanoparticles has made it possible to develop rapid immunoassays that offer far greater sensitivity and quantitative results.




The europium chelate complex possesses a longer fluorescence lifetime (μs) than traditional fluorophores (ns), allowing signal to be collected beyond the lifetime of background fluorescence. Moreover, its long Stokes shift ensures that incident light from the excitation source ($\lambda \sim 330\text{-}340\text{nm}$) does not interfere with collection of light by the detector ($\lambda \sim 610\text{-}620\text{nm}$). These traits, coupled with the availability of small, portable time-resolved fluorescence (TRF) readers, present new opportunities in the evolution of rapid diagnostics.

Our highly-fluorescent europium(III) nanoparticles offer exceptional functionality and stability for the development of diagnostic reagents. They have been utilized to develop highly sensitive assays in lateral flow and microplate formats, and are compatible with commercial europium chelate TRF readers. They are supplied in aqueous suspension at 1% solids (10mg/mL). See datasheet 741 for additional product details. The expiration date is 24 months from the date of shipment.

Comparison Data

As shown below, Bangs' europium chelate microspheres exhibit exceptional fluorescence intensity levels when compared to the leading competitor (1:60,000 dilutions, 350nm excitation, 610nm detection).



References

Bao DT, Kim DTH, Park H, Cuc BT, Ngok NM, Linh NTP, et al. (2017) *Rapid detection of avian influenza virus by fluorescent diagnostic assay using an epitope-derived peptide*. *Theranostics*; 7(7):1835-1846.

Liang R-L, Deng Q-T, Chen Z-H, Xu X-P, Zhou J-W, Liang J-Y, et al. (2017) *Europium (III) chelate microparticle-based lateral flow immunoassay strips for rapid and quantitative detection of antibody to hepatitis B core antigen*. *Sci Rep*; 7:14093.

Soukka T, Paukkunen J, Härmä H, Lönnberg S, Lindroos H, & Lövgren T. (2001). *Supersensitive time-resolved immunofluorometric assay of free prostate-specific antigen with nanoparticle label technology*. *Clin Chem*; 47(7), 1269-1278.

Tang Y, Zhang H, Liu X, Tromfimchuk E, Feng S, Ma T, et al. (2017) *Advantage of Eu³⁺-doped polystyrene microspheres compared with colloidal gold used in immunochromatographic assays for the detection of melamine in milk*. *J Food Sci*; 82(3):694-697.

Yeo S-J, Bao DT, Seo G-E, Bui CT, Kim DTH, Anh NTV, et al. (2017) *Improvement of a rapid diagnostic application of monoclonal antibodies against avian influenza H7 subtype virus using Europium nanoparticles*. *Sci Rep*; 7:7933.

EUROPIUM CHELATE

Cat. #	Product Description
FCEU001	0.10µm Europium Chelate COOH
FCEU002	0.20µm Europium Chelate COOH
FCEU003	0.30µm Europium Chelate COOH
21960	Europium Chelate COOH Sampler Pack 1mL of 0.10µm, 0.20µm, 0.30µm



Bangs Laboratories manufactures polymeric, silica and magnetic microsphere products setting the standards for diagnostic, research, and flow cytometry applications. No matter the project, we have a product that serves or we'll work to custom-design a solution to fit. And that's not the half of it.


We also stand behind our products. Regardless of the size of your question or the size of your company, we offer tech support, absolutely free.

Sound interesting? 

Visit: www.bangslabs.com

 [@particledoc](https://twitter.com/particledoc)

 info@bangslabs.com

 800.387.0672