



DESCRIPTION

Our new StarLight™ Calibration Slides feature vibrant ~6µm fluorescent microspheres dyed with a single fluorophore for basic imaging checks and calibration of fluorescence microscopes. Our four standard versions are appropriate for use with common microsphere filter sets: Glacial Blue (360, 450), Dragon Green (480, 520), Envy Green (525, 565), and Flash Red (660, 690), and are available individually or as a the full 4-slide collection.



Images of StarLight™ Calibration Slides. From Left to right : Glacial Blue, Dragon Green, Envy Green, and Flash Red.

SPECTRA



STORAGE AND STABILITY

Store at room temperature. Protect from light and moisture. StarLight™ slides do not have a formal expiration date. Longevity will depend on proper storage, as well as specific imaging conditions.

This product is for research use only and is not intended for use in humans or for *in vitro* diagnostic use.

ORDERING INFORMATION

| Cat. Code | Description | Sizes |
|-----------|---------------------------------------------|----------|
| SL1GB | StarLight™ Calibration Slide – Glacial Blue | 1 slide |
| SL1DG | StarLight™ Calibration Slide – Dragon Green | 1 slide |
| SL1EG | StarLight™ Calibration Slide – Envy Green | 1 slide |
| SL1FR | StarLight™ Calibration Slide – Flash Red | 1 slide |
| SL04K | StarLight™ Collection – Slide 4-Pack | 4 slides |

Order online anytime at www.bangslabs.com.